Q500 Sonicator

QSONICA

The Q500 is a powerful ultrasonic processor featuring programmable operation and a digital display of operating parameters. Popular applications include nanoparticle dispersion, creating emulsions, cell lysis and homogenization.

Adjustable pulse On and Off times can be programmed from 1 second to 1 minute. Total programming has a maximum setting of 10 hours. A wide variety of probes and accessories are available to handle virtually any application.





FEATURES:

Programmable operation Set time and amplitude for hands free operation

Pulse mode Prevent heat buildup in temperature sensitive samples

Digital amplitude / intensity control Output intensity can be set from 20-100%

Elapsed time indicator Displays duration of sonication

Display of wattage and joules Real-time energy monitoring

Overload protection Prevents damage to circuitry if a fault occurs

RoHS compliant

Uses lead free components

PART NO. Q500 INCLUDES:

- Generator
- Converter
- 1/2" diameter probe
- Power cable
- Converter cable
- Wrench set

TECHNICAL SPECIFICATIONS:

Power Rating:	500 watts
Frequency:	20 kHz
Programmable Timer:	10 hours
Adjustable Pulse On/Off:	1 second to 1 minute
Dimensions:	8"W x 15.25"L x 8.5"H
Voltage:	110V, 50/60Hz

Specify desired voltage for export.

Cup Horn

A Cup Horn offers indirect sonication and functions as a high intensity ultrasonic water bath. Multiple samples can be processed in sealed tubes eliminating cross contamination or aerosol issues.

The horn is mounted within an acrylic cup and the cup is filled with water. Sample tubes are placed in a rack at a fixed distance above the ultrasonic horn. Cavitation is produced in the water, processing the samples within the tubes. The #440 tube rack is included with the Cup Horn. This rack is made for 1.5mL polystyrene tubes which are proven to process samples more efficiently than 1.5mL polypropylene tubes.

Sonication generates heat so ports for cooling are located on each side of the cup. The #4905 Chiller is recommended for maintaining both the water temperature and water level within the Cup Horn.

The Sound Enclosure is highly recommended for all Cup Horn users. In addition to reducing sonication noise to safe levels, it securely holds the Cup Horn in place. Custom tube racks are available. Contact Qsonica for application assistance.



#431C2



Part #Tube Holder Description4408 x 1.5mL Polystyrene tubes4518 x 1.5mL Polypropylene tubes44912 x 600µl PCR tubes44524 x 300µl PCR tubes4551 x 50mL tube4541 x 15mL tube



QSONICA





Q700 / Q500 Accessories Flocells

Low Volume Flocells

The Low Volume Flocell (LVF) is available for use with either the Q500 or Q700 system. The Flocell (#4650) is equipped with ¼" (6 mm) hose barb fittings and does not include a probe. A ½" (12 mm) replaceable tip probe (#4643) or ½" (13 mm) solid probe (#4644) must be ordered separately. These probes feature a flange for proper mounting with the LVF. The replaceable tip probe is for use with aqueous samples only. Solid tip probes can be used with all types of solvents or low surface tension liquids.

Sonication generates heat so a Water Jacket (#4655) is available if the process requires cooling. The water jacket slides over the LVF and is used to recirculate cold water around the exterior of the flocell body. The water jacket includes 1/4" (6 mm) hose barb fittings.

The LVF is recommended for processing sample volumes above 1L. Routine applications include cell lysis, mixing, solubilizing and deagglomerating/ dispersing nanoparticles.







How Flocells Work

Flocells offer inline or continuous, large volume, batch sample processing. Flocells are ideal for mixing and dispersing applications. Batch volumes can be re-circulated through the system multiple times if increased sonication time is needed. Multiple units can be used in series to reduce processing time and/or maintain an even higher flow rate.

The liquid sample is pumped into the Flocell through the inlet at the bottom of the unit. As the sample passes through the cavitation field, it is processed. The processed liquid exits the unit through an outlet port. The degree of processing can be controlled by adjusting the intensity of sonication as well as flow rate.

High Volume Flocells

The High Volume Flocell (HVF) is available for use with either the Q700 or Q1375 System. The Flocell (#4549) is equipped with ½" (13 mm) sanitary connections, a water jacket and 1" (25 mm) Diameter probe (#4625). The water jacket can be used to recirculate cold water around the exterior of the flocell body. This helps reduce the heat generated during ultrasonic processing.

The HVF is recommended for processing batch volumes of 5L or more. Routine applications include cell lysis, mixing, solubilizing and deagglomerating/ dispersing nanoparticles.





Q700/Q500 Accessories Sound Enclosure



#432B2 (Probe sold separately.)

Part #	Description			
	Sound Enclosure with Converter Holder,			
432B2	Exterior Dimensions (W x H x D)			
	13.5 x 30.5 x 13 in.			
	(343 x 775 x 330 mm)			

Sound Enclosure

Sonicators are extremely loud devices and will cause discomfort to the user and anyone nearby. The Sound Enclosure reduces noise by approximately 20 dBa and is made to work with all accessories (excluding the Microplate Horn which has its own dedicated enclosure).

In addition to reducing noise, the Sound Enclosure has an internal support rod and converter mounting system. Any Qsonica probe or horn will be held safely and securely inside the unit.

Two ports are located on either side of the enclosure for coolant tubing or a temperature monitoring probe. The interior walls are lined with acoustical foam and the door has a window so experiments can be visually monitored.



(Cup horn sold separately.)



General Accessories

Replacement Converter





Temperature Monitoring Options



Replacement Converter Cable



K4

K4-10

6 ft. Long

10 ft. Long

Footswitch



Heavy Duty Stand

4130 Heavy Duty Stand The stand securely holds horns in place and ensures all tips are level. An adjustable jack stand to raise and lower samples is included Image: Clark stand stand with 1/2" Diameter Support Rod and Converter Clamp Part # Description 459 Stand with 1/2" Diameter Support Rod and Converter Clamp Jack Stand Stand with 1/2" Diameter Support Rod and Converter Clamp Jack Stand Baises and lower sample vessels a stationary prof as needed.	Part #	Descript	tion			
The stand securely holds horns in place and ensures all tips are level. An adjustable jack stand to raise and lower samples is included	4130	Heavy D Stanc	uty I			
Variable Description Variable Stand with V2" Diameter Support Rod Support Rod and converter Clamp Clamp	The stand securely holds horns in place and ensures all tips are level. An adjustable jack stand to raise and lower samples is included					
Large Clamp Stand Part # Description 459 Stand with ½" Diameter Support Rod and Converter Clamp Clamp Clamp Jack Stand Stand vith Kaises and lower Stationary profile Stationary profile Stationary profile Part # Description			2	5	2 27 27	
Part # Description 459 Stand with ½" Diameter Support Rod and Converter Clamp Jack Stand Raises and lowe sample vessels a stationary profase Part # Description	Large	Clamp	Sta	nd		
459 Stand with ½" Diameter Support Rod and Converter Clamp Jack Stand Raises and lower sample vessels ta a stationary prof as needed. Part # Description	Part #	Descrip	tion			
Jack Stand Raises and lower sample vessels a stationary prof as needed. Part # Description	459	Stand w 1⁄2" Diam Support and Conv Clam	vith leter Rod lerter p			-
Jack Stand Raises and lower sample vessels a stationary prof as needed. Part # Description						
Raises and lower sample vessels a stationary prof as needed. Part # Description	Jack \$	Stand				
Part # Description				Rais san a st as r	ses and aple ves ationary aeded.	lower sels to prob
			Part	#	Descr	iptio



High Throughput Horns

4 Tip Horns

The 4 Tip Horn enables 4 samples to be processed simultaneously. This horn is effective for cell disruption, mixing, homogenization and many other applications. This horn is available with either 1/8" or 1/4" tip diameters.

#4659 is a 4 Tip Horn with 1/8" tips. It is recommended for 1-15ml volumes and fits into both 1.5ml and 15ml tubes. The space between each tip is 1.05". The #504 and #510 Coolracks are compatible with #4659.

#4674 is a 4 Tip Horn with 1/4" tips. It is recommended for 10-50ml volumes and is

used with 15 and 50ml tubes. The space between each tip is 1.6". The #511 Coolrack is compatible with the #4674.



24 Tip Horn

The 24 Tip Horn processes each well of a 24 well plate at the same time. This horn is effective for cell disruption, mixing, dissolution and many other applications. The distance between each tip is .708".



#4597

Dual Horn

consistent results.

The #4130 Stand with Lab Jack is recommended for use with all high throughput horns. This stand holds the horn stable and level to ensure

The Dual Horn allows a single Sonicator unit to process two samples simultaneously. Two standard probes are attached to a rectangular shaped horn. The distance from center to center of each probe is 4.5". 3/4" solid tip probes are included with the Dual Horn but 1/2" or 1" probes can also be substituted.





Q700 Accessories

Indirect Horn Options

Microplate Horn

(Only for use with Q700)

Similar to a Cup Horn, but larger, the Microplate Horn is an indirect sonication device capable of processing an entire 96 well microtiter plate or many microtubes at one time.

Simply place your samples within the water-filled reservoir and the sonic energy is transferred into each individual well or tube.

The Horn is equipped with a clear acrylic collar to contain the liquid media within the reservoir. This allows the user to process deep well microplates or other tall vessels. Standard microtiter plates or PCR tubes require a smaller volume of liquid for sonication. For these applications, the clear acrylic collar may be removed and the lower, gray collar will allow for easier access to the samples.



Part #	Description		
Q700MPX	Q700 (Without Standard Probe), and the 431MPX		
431MPX	Microplate Horn, Pinch Clamps, Tubing and Sound Enclosure		
431MPXH	Microplate Horn Only		
432MP	Sound Enclosure for Microplate Horn		
444	300µl Microcentrifuge Tube Holder/Cover		

#Q700MPX

(with Chillerand tubing set)

GSONIC

Exterior dimensions of the Sound Enclosure are (W x D x H): 10 x 10 x 17 in. (254 x 245 x 432 mm).



The Microplate Horn is commonly used in PMCA research. A microcentrifuge tube holder and cover (#444) are available and often used for this application.



Q700 / Q500 Accessories

Direct Horn Options

Horns (also known as probes) are made from titanium and machined to specific sizes and shapes. When driven at their resonant frequency, they expand and contract longitudinally. This mechanical vibration is amplified and transmitted down the length of the probe. In liquid, the probe causes cavitation which constitutes the main mechanism for sample processing.

Choosing the appropriate horn is extremely important. The sample volume to be processed is directly related to the tip diameter. Smaller tip diameters (Microtip probes) deliver high intensity sonication, but the energy is focused within a small, concentrated area. Larger tip diameters can process larger volumes, but offer lower intensity. Probes are offered with replaceable, solid or sapphire tips.

Probe tips will pit or erode over time and require replacement. Replaceable tip probes are used with aqueous samples only. In addition to aqueous samples, Solid probes can be used with organic solvents, alcohols and low surface tension liquids. Contact Qsonica with questions regarding proper tip selection.



Standard Probes



Part #	Type of Tip	Processing Volume	Tip Diameter	Amplitude (microns)
4220	Replaceable Tip	20-250 ml	1/2" (13 mm)	120 µm
4219	Solid Tip	20-250 ml	1/2" (13 mm)	120 µm
4207	Replaceable Tip	50-500 ml	3/4" (19 mm)	60 µm
4208	Solid Tip	50-500 ml	3/4" (19 mm)	60 µm
4210	Replaceable Tip	100-1,000 ml	1" (25 mm)	30 µm
4209	Solid Tip	100-1,000 ml	1" (25 mm)	30 µm

info@sonicator.com • Tel. 1.877.338.9636

Q700/Q500 Accessories

Direct Horn Options

Replacement Tips for Standard Probes

Standard ½", ¾" and 1" horns have replaceable tips. During normal use, tips erode and become less effective over time. These worn tips can be easily removed and replaced.



Part #	Tip Diameter	For Use With
4406	1/2" (13 mm)	#4220
4407	3/4" (19 mm)	#4207
4408	1" (25 mm)	#4210





New Tip

Worn Tip

Microtip Probes

Microtips are thin, high intensity probes which are designed for processing small sample volumes. Microtips screw into the threaded end of the standard ½" probe (#4220).



Part #	Processing Volume*	Tip Diameter	Amplitude (microns)
4417	0.2-5 ml	1/16" (2 mm)	320 µm
4418	1-15 ml	1/8" (3 mm)	380 µm
4420	5-50 ml	1/4" (6 mm)	200 µm

* The recommended processing volume range is application specific. For example, samples containing surfactants foam easily and may require larger volumes for effective sonication. Contact Qsonica for application assistance.

Coupler with Stepped Microtip

The stepped microtip and coupler assembly is a low intensity option which can be used to process small volumes that do not require high power. The probe tip remains 1/8" in diameter for 48mm. This 2-piece assembly attaches directly to the converter.

Coupler — # 4421			Stepped Microtip — Pa	⊇ ırt # 4422
	Part #	Processing Volume	Tip Diameter	Amplitude (microns)
	4422*	0.5-15 ml	1/8″ (3 mm)	200 µm
	4421	Coupler - *requi	red for use of a Stepp	ed Microtip



Q700 / Q500 Accessories Direct Horn Options

Extenders

Standard probes may not be long enough to fit down into certain long necked vessels. Extender probes attach to standard horns of the same tip diameter and extend the length of the horn assembly. Extenders are available in 5" and 10" lengths with either solid, or replaceable tips.



Extenders offer the same processing volume and amplitude of their corresponding standard horn.

Part #	Type of Tip	Length	Tip Diameter
406HW	Solid Tip	5"	1/2" (13 mm)
406HWT	Replaceable Tip	5"	1/2" (13 mm)
407HW	Solid Tip	5"	3/4" (19 mm)
407HWT	Replaceable Tip	5"	3/4" (19 mm)
408HW	Solid Tip	5"	1" (25 mm)
408HWT	Replaceable Tip	5"	1" (25 mm)
407FW	Solid Tip	10"	3/4" (19 mm)
407FWT	Replaceable Tip	10"	3/4" (19 mm)
408FW	Solid Tip	10"	1" (25 mm)
408FWT	Replaceable Tip	10"	1" (25 mm)

Boosters



Booster horns increase the intensity of standard ³/₄" and 1" horns. Boosters attach between the converter and horn to increase amplitude by the gain ratio indicated below.

Part #	For Use With	Gain Ratio
4121	3/4" (19 mm) and 1" (25 mm) Probes	2 to 1

High Gain Horns



High gain horns (also known as high intensity horns) offer double the amplitude of standard ³/₄" and 1" horns. High gain horns attach directly to the converter.

Part #	Type of Tip	Processing Volume	Tip Diameter	Amplitude (microns)
4305	Replaceable Tip	50-500 ml	3/4" (19 mm)	120 µm
4306	Solid Tip	50-500 ml	3/4" (19 mm)	120 µm
4310	Solid Tip	100-1,000 ml	1" (25 mm)	60 µm
4311	Replaceable Tip	100-1,000 ml	1" (25 mm)	60 µm



Instrumentos Científicos

Anyover

Sonicadores

Reactores

Análisis Termico

Quimisorción

TPD/TPR/TPO

Porosidad

Espectrometría de Masas

Propiedades de Barrera

SIMS

Área Superficial BET

www.anyover.com.mx



En Anyover Instrumentación Científica distribuimos instrumentos científicos de alta tecnología para el campo de caracterización de materiales haciendo énfasis en la atención al cliente en temas de instalación, capacitación y soporte técnico. Nuestros representantes te apoyaran en la adquision y mantenimiento de instrumentos de nuestras marcas.



Thermo Fisher

Caracterización completa de materiales: Area Superficial BET, Porosidad, Quimisorción, Picnometría de Helio, TPD/TPR/TPO.





Espectrometros de masas para análisis de gases, sistemas de vacío, plasmas y ciencias de superficies.





Diseña, desarrolla y fabrica los sonicadores más avanzados del mundo.





Instrumentos para medir la permeabilidad al oxígeno (O2TR), al vapor de agua (WVTR) y al dióxido de carbono (CO2TR) y el innovador "data logger Wireless".





POROMETER

Fabricante belga alemán de equipos para caracterización de membranas y filtros.





Reactores para caracterización de catalizadores y otras aplicaciones a nivel laboratorio, para uso industrial y de investigación.





Instrumentos para Análisis Térmico (TGA, DSC, DTA, STA)

Sigue equipando tus investigaciones con nosotros.



Buscar



1



¿Cómo es estar con nosotros?



Comprar



3

Tengo el gusto de conocer al fundador de **ANYOVER** como un proveedor de confianza desde hace 10 años.

Mis más recientes adquisiciones: el Micro Reactor CATLAB y el Espectrómetro de Masas QGA, han sido de gran ayuda para nuestro grupo de investigación.



to Nos permiten flexibilidad en el diseño de nuestros experimentos así como confiabilidad en los resultados obtenidos y un incremento en el número de pruebas realizadas.

Dr. Heriberto Pfeiffer Instituto de materiales, UNAM ANYOVER Instrumentación Científica SA de CV Ciudad de México Tel. (55) 7676-5900 Contacto: ventas@anyover.com.mx www.anyover.com.mx



